

LIFE LINK TUBE FOR MOBI

Réf. 513720



LIFE LINK TUBE FOR MOBI



Réf. 513720

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe I
Marquage			Barcode 3700446028046
Guarantee (month):	12	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,199
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:50 - This article is available on our website : www.gsh-med.fr with the following reference 513720

All characteristics are given for information only and may be changed