

## 12 GLASSES FOR SPENGLER MANOMETERS VAQUEZ-LAUBRY CLASSIC/LIAN CLASSIC



Réf. 515273



# 12 GLASSES FOR SPENGLER MANOMETERS VAQUEZ-LAUBRY CLASSIC/LIAN CLASSIC



Réf. **515273**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446031787
Guarantee (month):	60	Customs Code	901890100000L
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	5,0
Length (cm) :		Length packaging (cm) :	5,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,052
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:28 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 515273

All characteristics are given for information only and may be changed