

SUNTECH CT50 BP STANDARD SPO2 TEMP PRINTER

Réf. DSU222204

SunTech Medical®

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0413

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE II



SUNTECH CT50 BP STANDARD SPO2 TEMP PRINTER

Réf. DSU222204

FEATURES				
Medical device	Yes	Classe IIa	Accuracy measurement	Classe II
Marquage	CE0413		Barcode	3661040214238
Guarantee (month):	24		Customs Code	
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	32,5
Length (cm) :			Length packaging (cm) :	37,5
Height (cm) :			Height packaging (cm) :	33,0
Weight (kg) :			Weight packaging (kg) :	5,644
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 02:32 - This article is available on our website : www.gsh-med.fr with the following reference DSU222204

All characteristics are given for information only and may be changed