

## DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM) CORAL

Réf. SPG445P



## DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM) CORAL



Réf. **SPG445P**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3700446041588
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	3,730
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:00 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SPG445P

All characteristics are given for information only and may be changed

## DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM) CORAL



**Réf. SPG445P**

### COMPOSITION

Code	Descriptive	Quantity
SPG445	AUTOTENSIO Upper arm (M/L) Blood Pressure Monitor Spengler Coral	6