

SET OF SILICONE TUBES Ø8X14 JAR + CONNECTOR FOR ASKIR C30BR

Réf. BCASPI115

CA-MI



CE
0123



SET OF SILICONE TUBES Ø8X14 JAR + CONNECTOR FOR ASKIR C30BR

CA-MI

Réf. BCASPI115

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	3661040151069
Guarantee (month):	12	Customs Code	90330090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,256
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:02 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI115

All characteristics are given for information only and may be changed

2/2