

RICHARDSON RETRACTOR 25 CM 44 X 38 MM

Réf. IE14439



RICHARDSON RETRACTOR 25 CM 44 X 38 MM



Réf. **IE14439**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040138183
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	10,8
Length (cm) :		Length packaging (cm) :	31,5
Height (cm) :		Height packaging (cm) :	6,0
Weight (kg) :	0,099	Weight packaging (kg) :	0,104
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:48 - This article is available on our website : www.gsh-med.fr with the following reference IE14439

All characteristics are given for information only and may be changed