

METAL TONGUE DEPRESSOR STAINLESS STEEL ADULT CURVED 135° 16 CM

Réf. IA02402



METAL TONGUE DEPRESSOR STAINLESS STEEL ADULT CURVED 135° 16 CM



Réf. IA02402

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040005966
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,2	Width packaging (cm):	
Length (cm) :	16,5	Length packaging (cm) :	
Height (cm) :	9,2	Height packaging (cm) :	
Weight (kg) :	0,066	Weight packaging (kg) :	0,066
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:40 - This article is available on our website : www.gsh-med.fr with the following reference IA02402

All characteristics are given for information only and may be changed