

## REAMER WITH "T" HANDLE SIZE 6 TO 14 MM (SET OF 9 PCS)

Réf. IA09300



## REAMER WITH "T" HANDLE SIZE 6 TO 14 MM (SET OF 9 PCS)



Réf. IA09300

FEATURES			
Medical device	No		
Marquage	Barcode	3661040080475	
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,450	Weight packaging (kg) :	0,450
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:44 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IA09300

All characteristics are given for information only and may be changed