

JACKSON TRACHEAL TUBE 3 PIECES N° 4 (DIA. 8 MM)

Réf. IC42304



JACKSON TRACHEAL TUBE 3 PIECES N° 4 (DIA. 8 MM)



Réf. IC42304

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040011585
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,5	Width packaging (cm):	
Length (cm) :	8,4	Length packaging (cm) :	
Height (cm) :	2,0	Height packaging (cm) :	
Weight (kg) :	0,030	Weight packaging (kg) :	0,030
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:14 - This article is available on our website : www.gsh-med.fr with the following reference IC42304

All characteristics are given for information only and may be changed