

ELECTRODES FOR TRAINING DEFIBRILLATOR (5 PAIRS)

Réf. SJACCE69

SaverOne®



ELECTRODES FOR TRAINING DEFIBRILLATOR (5 PAIRS)



Réf. **SJACCE69**

FEATURES			
Medical device	No		
Marquage		Barcode	3661040137285
Guarantee (month):	12	Customs Code	90189085
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	16,0	Width packaging (cm):	
Length (cm) :	23,0	Length packaging (cm) :	
Height (cm) :	22,0	Height packaging (cm) :	
Weight (kg) :	0,100	Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:37 - This article is available on our website : www.gsh-med.fr with the following reference SJACCE69

All characteristics are given for information only and may be changed