

LID FOR 5000ML POLYCARBONATE JAR NEW HOSPIVAC

Réf. BCASPI120

CA-MI



LID FOR 5000ML POLYCARBONATE JAR NEW HOSPIVAC



Réf. BCASPI120

FEATURES			
Medical device	No		
Marquage		Barcode	000000000710
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	18,5
Length (cm) :		Length packaging (cm) :	18,5
Height (cm) :		Height packaging (cm) :	12,0
Weight (kg) :		Weight packaging (kg) :	0,272
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:35 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI120

2/2