

## AUTOCLAVABLE SILICONE CUFF, 61 X 11CM (24 INCH), SINGLE

Réf. 554072



## AUTOCLAVABLE SILICONE CUFF, 61 X 11CM (24 INCH), SINGLE



Réf. **554072**

FEATURES			
Medical device	No		
Marquage	CE	Barcode	28809315676571
Guarantee (month):	12	Customs Code	9018 90 2090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	11,0	Width packaging (cm):	
Length (cm) :	61,0	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:50 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 554072

All characteristics are given for information only and may be changed