



## CONVERTIBLE BACK PACK ETRETAT, DE BOISSY, WITH BROWN LEATHER POCKET

*De Boissy*  
— DEPUIS 1982

Réf. MDETRE03



## CONVERTIBLE BACK PACK ETRETAT, DE BOISSY, WITH BROWN LEATHER POCKET



Réf. **MDETRE03**

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040218984
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	15,0	Width packaging (cm):	50,0
Length (cm) :	30,0	Length packaging (cm) :	50,0
Height (cm) :	45,0	Height packaging (cm) :	35,0
Weight (kg) :	1,860	Weight packaging (kg) :	1,860
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



03/05/2022 21:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference MDETRE03

All characteristics are given for information only and may be changed