

KOERTE RETRACTOR 8 PRONGS 26X40MM BLUNT/250MM

Réf. IE11501



KOERTE RETRACTOR 8 PRONGS 26X40MM BLUNT/250MM



Réf. IE11501

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040044187
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,0	Width packaging (cm):	
Length (cm) :	24,6	Length packaging (cm) :	
Height (cm) :	3,0	Height packaging (cm) :	
Weight (kg) :	0,107	Weight packaging (kg) :	0,107
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:10 - This article is available on our website : www.gsh-med.fr with the following reference IE11501

All characteristics are given for information only and may be changed