

KILLIAN NASAL SPECULUM 30 MM WITH LATERAL SCREW

Réf. IS23030



KILLIAN NASAL SPECULUM 30 MM WITH LATERAL SCREW



Réf. IS23030

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040031149
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,0	Width packaging (cm):	
Length (cm) :	14,8	Length packaging (cm) :	
Height (cm) :	3,0	Height packaging (cm) :	
Weight (kg) :	0,101	Weight packaging (kg) :	0,101
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:46 - This article is available on our website : www.gsh-med.fr with the following reference IS23030

All characteristics are given for information only and may be changed