

KILLIAN NASAL SPECULUM 50 MM WITH LATERAL SCREW

Réf. IS23050



KILLIAN NASAL SPECULUM 50 MM WITH LATERAL SCREW



Réf. **IS23050**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040031156
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,5	Width packaging (cm):	
Length (cm) :	15,0	Length packaging (cm) :	
Height (cm) :	5,2	Height packaging (cm) :	
Weight (kg) :	0,099	Weight packaging (kg) :	0,099
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:46 - This article is available on our website : www.gsh-med.fr with the following reference IS23050

All characteristics are given for information only and may be changed