

LEBSCHER STERNUM CHISEL 24 CM

Réf. IS24424



LEBSCHE STERNUM CHISEL 24 CM



Réf. IS24424

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040031224
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,5	Width packaging (cm):	10,5
Length (cm) :	24,0	Length packaging (cm) :	24,0
Height (cm) :	2,0	Height packaging (cm) :	2,0
Weight (kg) :	0,214	Weight packaging (kg) :	0,214
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:17 - This article is available on our website : www.gsh-med.fr with the following reference IS24424

All characteristics are given for information only and may be changed