

## JARS POLYCARBONATE 4L WITH LID FOR HOSPIVAC 350/400

CA-MI

Réf. BCASPI82

## JARS POLYCARBONATE 4L WITH LID FOR HOSPIVAC 350/400

CA-MI

Réf. BCASPI82

FEATURES			
Medical device	No		
Marquage		Barcode	3661040169323
Guarantee (month):	0	Customs Code	90330090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:33 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI82