

## DISSECTING SET

Réf. IT01600



# DISSECTING SET



Réf. **IT01600**

Consisting of:

- IP26410 Splinter forceps, 10 cm
- IP26014 Dissection forceps without teeth, 14 cm
- IS02414, Fluted probe, 14 cm
- ID07601, Dissecting needle with handle, 14 cm
- ID08601, Dissecting cutting needle with handle, 14 cm
- IC11414 Straight round-ended scissors, 14 cm
- IC09411 Straight iridectomy scissors, 11.5 cm
- ID05404, no.4 scalpel handle

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040031569
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	13,5	Width packaging (cm):	25,0
Length (cm) :	18,0	Length packaging (cm) :	15,0
Height (cm) :	3,0	Height packaging (cm) :	
Weight (kg) :	0,282	Weight packaging (kg) :	0,282
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:48 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IT01600

All characteristics are given for information only and may be changed