

## INSTRUMENT TROLLEY S.S 3 SHELVES WITH SIDE RAILS 60X40X80 CM

Réf. OSA6041



# INSTRUMENT TROLLEY S.S 3 SHELVES WITH SIDE RAILS 60X40X80 CM



Réf. OSA6041

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040168296
Guarantee (month):	12	Customs Code	9402900000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	40,0	Width packaging (cm):	45,5
Length (cm) :	60,0	Length packaging (cm) :	85,0
Height (cm) :	80,0	Height packaging (cm) :	14,0
Weight (kg) :	12,200	Weight packaging (kg) :	12,200
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/10/2022 12:38 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference OSA6041

All characteristics are given for information only and may be changed