

ALUMINIUM ALLOY FOLDAWAY DISASTER STRETCHER, BLUE

Réf. SBRAN07_B



ALUMINIUM ALLOY FOLDAWAY DISASTER STRETCHER, BLUE



Réf. SBRAN07_B

- Stretcher foldable in 2, with wheels and stabilisers
- High resistance aluminium structure and waterproof fabric
- The back of this stretcher can be raised up to 60°
- Lightweight, easy to carry
- Dimensions (L x W x H): 185 X 50 x 24 cm, unfolded; 92 x 50 x 11 cm, folded
- Weight: 8 kg
- Maximum load carried: 159 kg
- Single colour: blue

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040094625
Guarantee (month):	12	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	50,0	Width packaging (cm):	12,0
Length (cm) :	185,0	Length packaging (cm) :	54,0
Height (cm) :	24,0	Height packaging (cm) :	98,0
Weight (kg) :	7,000	Weight packaging (kg) :	9,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:20 - This article is available on our website : www.gsh-med.fr with the following reference SBRAN07_B