

HEAD IMMOBILIZER, RED

Réf. SBRAN10_O

 Holtex

GARANTIE
1
AN
TIAS
GARANTIE

CE

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
I



HEAD IMMOBILIZER, RED



Réf. SBRAN10_O

- Manufactured with high density plastic materials, wide holes at the ears to monitor the auditory canal of the patient
- It can be used with a stretcher and a hard board
- Minimum Interference with radiography systems or MRI.
- Plastic Coated, water-resistant and easy to clean
- Dimensions (L x W x H): 40 x 25 x 16 cm
- Weight: 1 kg
- Single colour: orange

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040094656
Guarantee (month):	12	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	25,0	Width packaging (cm):	26,0
Length (cm) :	40,0	Length packaging (cm) :	43,0
Height (cm) :	16,0	Height packaging (cm) :	14,0
Weight (kg) :	1,000	Weight packaging (kg) :	0,950
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:21 - This article is available on our website : www.gsh-med.fr with the following reference SBRAN10_O