

## BECKMANN RETRACTOR (4X4 TEETH) 21 CM

Réf. IE13421



## BECKMANN RETRACTOR (4X4 TEETH) 21 CM



Réf. IE13421

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016634
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,0	Width packaging (cm):	
Length (cm) :	20,4	Length packaging (cm) :	
Height (cm) :	2,5	Height packaging (cm) :	
Weight (kg) :	0,123	Weight packaging (kg) :	0,123
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:11 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IE13421

All characteristics are given for information only and may be changed