

BECKMANN RETRACTOR 33 CM 3X4 TEETH 6 CM

Réf. IE13434



BECKMANN RETRACTOR 33 CM 3X4 TEETH 6 CM



Réf. IE13434

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016665
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,5	Width packaging (cm):	
Length (cm) :	33,5	Length packaging (cm) :	
Height (cm) :	6,0	Height packaging (cm) :	
Weight (kg) :	0,356	Weight packaging (kg) :	0,356
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:33 - This article is available on our website : www.gsh-med.fr with the following reference IE13434

All characteristics are given for information only and may be changed