

## NORMAL DELIVERY AND EPISIOTOMY REPAIR SET

Réf. IBOI117



## NORMAL DELIVERY AND EPISIOTOMY REPAIR SET



Réf. IBOI117

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe I                |               |
| Marquage                   | CE        | Barcode                 | 3661040081748 |
| Guarantee (month):         | 24        | Customs Code            | 9018.9090     |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   |               |
| Length (cm) :              |           | Length packaging (cm) : |               |
| Height (cm) :              |           | Height packaging (cm) : |               |
| Weight (kg) :              | 1,400     | Weight packaging (kg) : | 1,400         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



30/03/2022 15:15 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IBOI117

All characteristics are given for information only and may be changed

# NORMAL DELIVERY AND EPISIOTOMY REPAIR SET

Réf. IBO117

## COMPOSITION

| Code    | Descriptive                                  | Quantity |
|---------|--|----------|
| IA06520 | Reverdin needle - average curved 20 cm       | 1        |
| IC07418 | Mayo Scissors Straight 18 cm                 | 1        |
| IC07518 | Mayo Scissors Curved 18 cm                   | 1        |
| IC30108 | Bar Umbilical Clamp 8 cm                     | 1        |
| II05200 | Ribemont Insufflator (without bulb), 9cm     | 1        |
| AI25105 | S/S Box 25 x 10 x 5 cm                       | 1        |
| IP02016 | Kocher Forceps 1x2 Teeth Straight 16 cm      | 2        |
| IP02316 | Rochester-Pean Forceps Curved no teeth 16 cm | 1        |
| IP10614 | Backhaus Towel Forceps 14 cm                 | 2        |
| IP13624 | Cheron forceps 24 cm                         | 1        |
| IP15022 | J.L Faure Forceps Curved 1x2 Teeth 22 cm     | 1        |
| IP26018 | Tissue Forceps 1x2 Teeth 18 cm               | 1        |
| IP26216 | Dissecting Forceps No Teeth 16 cm            | 1        |
| IP44401 | Membrane Punch 19 cm                         | 1        |
| IS02414 | Director and Tongue tie 14 cm                | 1        |
| IS03402 | Metal Female Catheter 15 cm                  | 1        |