

CHILD SENSOR "DEPORTED" VELCRO FOR OXIMETER PC-60E

Réf. DPOXY05



CHILD SENSOR "DEPORTED" VELCRO FOR OXIMETER PC-60E



Réf. DPOXY05

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040165998
Guarantee (month):	6	Customs Code	9018909911
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	15,0
Length (cm) :		Length packaging (cm) :	14,5
Height (cm) :		Height packaging (cm) :	1,0
Weight (kg) :	0,059	Weight packaging (kg) :	0,059
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:12 - This article is available on our website : www.gsh-med.fr with the following reference DPOXY05

All characteristics are given for information only and may be changed