

## EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Réf. MEEMER02\_B

GARANTIE  
**2**  
ANS  
YEARS  
GARANTEE



## EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

**EMERGENCY'S™**

Réf. MEEMER02\_B

FEATURES			
Medical device	No		
Marquage	Barcode	8436564540874	
Guarantee (month):	24	Customs Code	420212
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	25,0	Width packaging (cm):	24,0
Length (cm) :	44,0	Length packaging (cm) :	30,2
Height (cm) :	27,0	Height packaging (cm) :	22,0
Weight (kg) :	1,290	Weight packaging (kg) :	1,300
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



12/04/2022 21:09 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference MEEMER02\_B

2/2