

EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLACK

Réf. MEEMER02_N

EMERGENCY'S™

GARANTIE
2
ANS
YEARS
GARANTEE



EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLACK

EMERGENCY'S™

Réf. MEEMER02_N

| FEATURES | | | |
|----------------------------|-----------|-------------------------|--------|
| Medical device | No | | |
| Marquage | Barcode | 8436564542977 | |
| Guarantee (month): | 24 | Customs Code | 420212 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 26,0 | Width packaging (cm): | 24,0 |
| Length (cm) : | 31,6 | Length packaging (cm) : | 30,2 |
| Height (cm) : | 19,0 | Height packaging (cm) : | 22,0 |
| Weight (kg) : | 1,290 | Weight packaging (kg) : | 1,300 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



12/04/2022 21:10 - This article is available on our website : www.gsh-med.fr with the following reference MEEMER02_N

2/2