

## CLINIC CUFF/1 TUBE/WITH COUPLING # 316 151 AND CLAMP/INFANT (XS)

Réf. 337846



## CLINIC CUFF/1 TUBE/WITH COUPLING # 316 151 AND CLAMP/INFANT (XS)



Réf. **337846**

Single tube Clinic® cuff for Lian® + clamp compatible w/ Spengler hand-held BPMs Protocol patient monitors

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446000134
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	10,0
Length (cm) :		Length packaging (cm) :	13,0
Height (cm) :		Height packaging (cm) :	4,5
Weight (kg) :		Weight packaging (kg) :	0,088
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:33 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 337846

All characteristics are given for information only and may be changed