

## PACK OF 3 CLINIC CUFFS/1 TUBE/ WITH CONNECTOR (L, M, S)

Réf. 337849



# PACK OF 3 CLINIC CUFFS/1 TUBE/ WITH CONNECTOR (L, M, S)



**Réf. 337849**

Single tube Clinic® cuffs for Lian® + clamp compatible w/ Spengler hand-held BPMs Protocol patient monitors

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe I                |               |
| Marquage                   |           | Barcode                 | 3700446030865 |
| Guarantee (month):         | 6         | Customs Code            | 9018901000    |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   | 17,0          |
| Length (cm) :              |           | Length packaging (cm) : | 24,0          |
| Height (cm) :              |           | Height packaging (cm) : | 8,0           |
| Weight (kg) :              |           | Weight packaging (kg) : | 0,250         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



26/04/2022 12:33 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 337849

All characteristics are given for information only and may be changed