

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - YELLOW (EM13.002)

EMERGENCY'S™

Réf. MEEMER02_J

GARANTIE
2
ANS
YEARS
GARANTEE



EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - YELLOW (EM13.002)

EMERGENCY'S™

Réf. **MEEMER02_J**

Light and practical emergency bag with two internal adjustable separators and a zipper pocket in the lid to facilitate the organisation of health equipment. It also has three external pockets; the front one is divided with elastic rubber bands.

- High resistance zip closures
- Reflective strips
- Washable briefcase
- Reinforced handles and straps
- Dimensions 44 x 25 x 27 cm; capacity 22 l
- Material 600D Polyester
- Weight 1.29 kg

FEATURES			
Medical device	No		
Marquage	Barcode	8436564540782	
Guarantee (month):	24	Customs Code	420212
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	26,0	Width packaging (cm):	26,0
Length (cm) :	31,6	Length packaging (cm) :	31,6
Height (cm) :	19,0	Height packaging (cm) :	19,0
Weight (kg) :		Weight packaging (kg) :	1,300
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



12/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference MEEMER02_J