

EMERGENCY LIGHT BAG, RED

Réf. MEEMER02_R

EMERGENCY'S™



GARANTIE
2
ANS
YEARS
GUARANTEE



EMERGENCY LIGHT BAG, RED

EMERGENCY'S™

Réf. **MEEMER02_R**

Light and practical emergency bag with two internal adjustable separators and a zipper pocket in the lid to facilitate the organisation of health equipment. It also has three external pockets; the front one is divided with elastic rubber bands.

- High resistance zip closures
- Reflective strips
- Washable briefcase
- Reinforced handles and straps
- Dimensions 44 x 25 x 27 cm; capacity 22 l
- Material 600D Polyester
- Weight 1.29 kg

FEATURES			
Medical device	No		
Marquage		Barcode	8436564540775
Guarantee (month):	24	Customs Code	420212
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	25,0	Width packaging (cm):	24,0
Length (cm) :	44,0	Length packaging (cm) :	30,2
Height (cm) :	27,0	Height packaging (cm) :	22,0
Weight (kg) :	1,290	Weight packaging (kg) :	1,300
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



12/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference MEEMER02_R

2 / 2

All characteristics are given for information only and may be changed