

EMERGENCY LIGHT BAG, RED

EMERGENCY'S™

Réf. MEEMER03_R

GARANTIE
2
ANS
YEARS
GARANTIE



EMERGENCY LIGHT BAG, RED

EMERGENCY'S™

Réf. **MEEMER03_R**

This bag is designed to carry the most fragile emergency equipment.

It has four compartments for internal supports and high capacity compartment.

Transparent pocket in the inside lid, designed in order to be able to see and to organise the material easily, thanks to the rubber bands which ensure correct attachment.

It also has five external pockets, two side pockets, one on the back and one on the top and a front pocket which can be divided internally thanks to pockets and rubber bands.

- Waterproof bottom bag
- High resistance zips
- Reflective strips
- Reinforced handles and straps
- External card holder
- Transparent top pocket
- Dimensions 55 x 31.5 x 33 cm; capacity 38 l
- Material 600D Polyester
- Weight 1.76 kg

FEATURES			
Medical device	No		
Marquage	Barcode	8436564540799	
Guarantee (month):	24	Customs Code	420212
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	31,5	Width packaging (cm):	28,0
Length (cm) :	55,0	Length packaging (cm) :	50,0
Height (cm) :	33,0	Height packaging (cm) :	20,0
Weight (kg) :	1,760	Weight packaging (kg) :	1,800
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



12/04/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference MEEMER03_R

2 / 2

All characteristics are given for information only and may be changed