

## LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)

Réf. 519001



# LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



**Réf. 519001**

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

| FEATURES                   |           |                                  |               |
|----------------------------|-----------|----------------------------------|---------------|
| Medical device             | Yes       | Classe Im (fonction de mesurage) |               |
| Marquage                   | CE0459    | Barcode                          | 3700446007225 |
| Guarantee (month):         | 60        | Customs Code                     | 90189010      |
| Sterile                    | No        |                                  |               |
| INFORMATIONS               | Main unit | Packaged unit                    |               |
| Width (cm) :               |           | Width packaging (cm):            |               |
| Length (cm) :              |           | Length packaging (cm) :          |               |
| Height (cm) :              |           | Height packaging (cm) :          |               |
| Weight (kg) :              |           | Weight packaging (kg) :          | 0,905         |
| SPECIFIC STORAGE CONDITION |           |                                  |               |
| Temperature min            | -20       | Temperature max                  | 70            |
| Humidity min               | 85        | Humidity max                     | 85            |



21/11/2022 21:09 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 519001

All characteristics are given for information only and may be changed