

## PULSE STETHOSCOPE / PEARL BLUE



Réf. 507524



CE



# PULSE STETHOSCOPE / PEARL BLUE



**Réf. 507524**

Pearly tube model. For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaural with built-in spring.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446003968
Guarantee (month):	36	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm) :	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,220
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



02/12/2022 12:29 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 507524

All characteristics are given for information only and may be changed