

LIAN CLASSIC / GREY COTTON CUFF / ADULT (M)

Réf. 519011



LIAN CLASSIC / GREY COTTON CUFF / ADULT (M)



Réf. 519011

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

| FEATURES | | | |
|----------------------------|-----------|----------------------------------|---------------|
| Medical device | Yes | Classe Im (fonction de mesurage) | |
| Marquage | CE0459 | Barcode | 3700446007218 |
| Guarantee (month): | 60 | Customs Code | 90189010 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | |
| Length (cm) : | | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | | Weight packaging (kg) : | 0,706 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | -20 | Temperature max | 70 |
| Humidity min | 85 | Humidity max | 85 |



26/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference 519011

All characteristics are given for information only and may be changed