

LIAN CLASSIC / CARTOON COTTON CUFF / CHILD (S)

Réf. 519029



LIAN CLASSIC / CARTOON COTTON CUFF / CHILD (S)



Réf. 519029

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446008734
Guarantee (month):	60	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,663
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	85	Humidity max	85



26/04/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference 519029

All characteristics are given for information only and may be changed