

## LIAN CLASSIC /MULTI COTTON CUFF/ (L, M, S)

Réf. 519050



CE  
0459



## LIAN CLASSIC /MULTI COTTON CUFF/ (L, M, S)



**Réf. 519050**

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446024789
Guarantee (month):	60	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	18,5
Length (cm) :		Length packaging (cm) :	26,5
Height (cm) :		Height packaging (cm) :	8,9
Weight (kg) :		Weight packaging (kg) :	1,150
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	85	Humidity max	85



26/04/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 519050

All characteristics are given for information only and may be changed