

## LIAN METAL / CARTOON DESIGN COTTON CUFF / NEWBORN (XXS)



Réf. 521314



CE  
0459



## LIAN METAL / CARTOON DESIGN COTTON CUFF / NEWBORN (XXS)



Réf. **521314**

Large dial, ambidextrous model (simply rotate to change hands). Metal casing with large dial for easy reading from all angles. Soft PVC manometer cover. Ergonomic spoon handle and inflation bulb. Cuffs with screw-on metal quick coupling. Supplied in a nylon pouch.

| FEATURES                   |           |                                  |               |
|----------------------------|-----------|----------------------------------|---------------|
| Medical device             | Yes       | Classe Im (fonction de mesurage) |               |
| Marquage                   | CE0459    | Barcode                          | 3700446006655 |
| Guarantee (month):         | 36        | Customs Code                     | 90189010      |
| Sterile                    | No        |                                  |               |
| INFORMATIONS               | Main unit | Packaged unit                    |               |
| Width (cm) :               |           | Width packaging (cm):            | 13,0          |
| Length (cm) :              |           | Length packaging (cm) :          | 19,0          |
| Height (cm) :              |           | Height packaging (cm) :          | 7,0           |
| Weight (kg) :              |           | Weight packaging (kg) :          | 0,520         |
| SPECIFIC STORAGE CONDITION |           |                                  |               |
| Temperature min            | -20       | Temperature max                  | 70            |
| Humidity min               | 85        | Humidity max                     | 85            |



26/04/2022 21:09 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 521314

All characteristics are given for information only and may be changed