

## LIAN METAL / MULTI CUFFS (L, M, S)



Réf. 521330



CE  
0459



## LIAN METAL / MULTI CUFFS (L, M, S)



**Réf. 521330**

Large dial, ambidextrous model (simply rotate to change hands). Metal casing with large dial for easy reading from all angles. Soft PVC manometer cover. Ergonomic spoon handle and inflation bulb. Cuffs with screw-on metal quick coupling. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446009939
Guarantee (month):	36	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	18,5
Length (cm) :		Length packaging (cm) :	26,5
Height (cm) :		Height packaging (cm) :	9,0
Weight (kg) :		Weight packaging (kg) :	1,213
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	85	Humidity max	85



26/04/2022 21:09 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 521330

All characteristics are given for information only and may be changed