

CHILD CUFF 2 TUBES BLACK FOR MOBILE



Réf. DS99123



CHILD CUFF 2 TUBES BLACK FOR MOBILE



Réf. **DS99123**

FEATURES			
Medical device	No		
Marquage	Barcode	3661040074825	
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,0	Width packaging (cm):	
Length (cm) :	15,1	Length packaging (cm) :	
Height (cm) :	2,0	Height packaging (cm) :	
Weight (kg) :	0,091	Weight packaging (kg) :	0,091
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:58 - This article is available on our website : www.gsh-med.fr with the following reference DS99123