

LARGE ADULT CUFF 2 TUBES BLACK FOR MOBILE

Réf. DS99125



CE
0483



LARGE ADULT CUFF 2 TUBES BLACK FOR MOBILE



Réf. **DS99125**

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0483	Barcode	3661040074849
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,450
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:59 - This article is available on our website : www.gsh-med.fr with the following reference DS99125

All characteristics are given for information only and may be changed