

10'3M PATIENT BP HOSE BAYONET TO BAYONET

Réf. DSU551456

 SunTech Medical®



10'3M PATIENT BP HOSE BAYONET TO BAYONET

Réf. DSU551456

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE		Barcode	3700446032562
Guarantee (month):	12		Customs Code	
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	2,0
Length (cm) :			Length packaging (cm) :	19,0
Height (cm) :			Height packaging (cm) :	18,0
Weight (kg) :			Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



05/04/2022 21:06 - This article is available on our website : www.gsh-med.fr with the following reference DSU551456