

## LIAN NANO / SILVER RING / BLACK NYLON / ADULT (M)

Réf. 513210



CE  
0459



# LIAN NANO / SILVER RING / BLACK NYLON / ADULT (M)



Réf. 513210

**Sphygmomanometer (adult), aneroid sphygmomanometer. Measures blood pressure non-invasively by displaying the pressure in a cuff wrapped around a patient's arm. The systolic and diastolic pressure is usually assessed by listening to Korotkoff sounds generated by arterial blood flow using a stethoscope simultaneously.**

#### HIGH PRECISION MEASUREMENTS

- High precision "watchmaking" mechanism

#### COMFORT

- Lightweight
- Ergonomic spoon and grooved bulb for perfect handling
- Large dial for easy reading
- Large inflation bulb for rapid inflation
- SPENGLER Velcro fastening, the perfect combination of maintenance and practicality

#### SHOCK RESISTANT

- Shock resistant ABS casing

#### CUFF

- Adult cuff size M (26-34 cm), Oxford Nylon, black

#### OPTION

- 6 cuff sizes (XXS, XS, S, M, L, XL) available and Clinic version to prevent cross-contamination

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446009953
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	13,5
Length (cm) :		Length packaging (cm) :	19,0
Height (cm) :		Height packaging (cm) :	7,0
Weight (kg) :	0,400	Weight packaging (kg) :	0,500
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



26/04/2022 21:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 513210

All characteristics are given for information only and may be changed