

LIAN NANO/SILVER RING/ CARTOON DESIGN COTTON / NEWBORN (XXS)

Réf. 513249



CE
0459



LIAN NANO/SILVER RING/ CARTOON DESIGN COTTON / NEWBORN (XXS)



Réf. **513249**

Lightweight and robust hand-held, combines style and performance. Lightweight, shock-resistant ABS casing. Easy-read large dial. Cartoon-pattern cotton cuffs with screw-on metal quick coupling. Ergonomic inflation bulb and brass valve with ABS cap. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446024307
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,478
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



26/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference 513249

All characteristics are given for information only and may be changed