

CB48

HEATING BLANKETS FOR OPERATING AND RECOVERY ROOM

TECHNICAL SHEET



MADE IN
FRANCE



- Reusable!
- Easy to use
- Safe
- Made In France

Clinical studies demonstrated that perioperative hypothermia may cause complications for patients. CB48 heating blankets are patient active warming systems dedicated to operative rooms, post-interventional monitoring rooms (SSPI) and for warming up patients prior to surgery.

CE
0459

GSH

30 rue Jean de Guiramand
13290 Aix en Provence | FRANCE
www.gsh-med.fr

Contact :
+33 (0) 4 42 90 31 31 - info@gsh-med.fr

CONTROL BOX TECHNICAL SPECIFICATIONS

WEIGHT		1.7 Kg
TEMPERATURE CONTROL	Measurement principle	CTN
	Measurement range	20°C to 41°C
	Measurement accuracy	+/- 0.8°C
	Adjustment range	30°C to 40°C / 0.5°C step
CONDITIONS OF USE	Recommended room temperature	21°C to 26°C
	Temperature range	10°C / 40°C
	Atmospheric pressure	700 hPa to 1060 hPa
	Hygrometry (non condensing)	30% to 90%
STORAGE CONDITIONS	Type of packaging	Origin
	Local	Ventilated
	Storage temperature	-20°C to 50°C
	Atmospheric pressure	900 hPa to 1060 hPa
PROTECTION	Hygrometry (non condensing)	30% to 90%
	BF Type / IP32	
CE MARKING	CE Marking 0459, notified body : G-MED, Device Class : IIb	

HEATING BLANKETS TECHNICAL SPECIFICATIONS

	CB48 HEATING BLANKET REFERENCE:		
	CB48-2 Model Pediatric- Half body	CB48-3 Model Adult - Upper body	CB48-11 Model Adult - Full Body
DIMENSIONS	80 x 100 cm	50 X 150 cm	100 x 150 cm
WEIGHT (w/out box) +/-0.2 Kg	1,9 Kg	1.7 Kg	2.7 Kg
POWER	70W	55 W	90 W
POWER SUPPLY	230 V-50 Hz		
HEATING VOLTAGE	48Vcc		
PROTECTION	BF		
WATERPROOFING	IP 35		
COVER	PU		
INNER FABRIC	83% acrylic / 17% cotton		
HANDLE	0	0	4
TEMPERATURE SETTING	30 to 40°C		
SKIN PROBE	no		

References

DESCRIPTION	REFERENCE
Heating Blanket - Half body/Pediatric	CB48-2+BCB48
Heating Blanket - Adult (half body)	CB48-3+BCB48
Heating Blanket - Adult (full body)	CB48-11+BCB48



Accessories:

DESCRIPTION	REFERENCE
Cover holder for serum stand	SUPPS



ECM FRENCH QUALITY

CE MARKING 0459

Warming cots are class IIb medical devices. They comply with current European regulations.

CUSTOMER TECHNICAL SUPPORT CENTRE

GSH provides maintenance and repair services for heating blankets. Contact our after-sales service on : +33 5 55 56 32 40.

ECM WARRANTY

ECM heating blankets are guaranteed for 12 months (under normal conditions of use). The warranty does not cover maintenance faults made by the user.