

CONE SINGLE 90X 12CM (35.4 INCH) NO PINCH CUFF (THIGH - L)



Ref. 554057



25/08/2024 18:35:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CONE SINGLE 90X 12CM (35.4 INCH) NO PINCH CUFF (THIGH - L)



Ref. 554057

Description

Single conical cuffs can be used for a wide spectrum of surgical procedures.

Armbands with zero-pinch technology: the innovative design of the withers of the withers DTS 3000 allows uniform and homogeneous inflation over the entire circumference of the cuff. Better wrapped, the limb requires less pressure to stop blood circulation, reducing here too the risk of complications.

Features

CE Marking	CE	Medical device	Class I
Bar code	28809315676250	HS Code	9018 90 2090
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	12
Packed unit / Height (cm)	nc	Packed unit / Length (cm)	nc
Bare Unit / Length (cm)	90	Packed unit / Weight	0.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	60
Minimum humidity	5	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 18:35:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr