

## DISPOSABLE FLOVAC LINER 2L FOR ASPIRET, ASKIR 20, 30, C30BR



Ref. BCASPI100



25/08/2024 18:27:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DISPOSABLE FLOVAC LINER 2L FOR ASPIRET, ASKIR 20, 30, C30BR



Ref. BCASPI100

### Description

- Flexible disposable polyethylene pouch

- Hermetically welded to a rigid lid inside which a hydrophobic and antibacterial anti-reflux filter is placed. This filter also works as an overflow, suction valve deactivation system whenever the maximum filling capacity is reached.

- The pocket must be placed in its FLOVAC® rigid polycarbonate autoclavable CONTAINER.

FLOVAC® systems therefore replace the use of autoclavable jars, which reduces cleaning and sterilization costs.

- Compatible with New Aspiret, Series Askir and Hospivac series Suction Units.

### Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3661040151045	HS Code	90189084
Warranty (months)	12	Country	Italy

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	14
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	5.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	17.5
Packed unit / Weight	0.2	Bare Unit / Weight (kg)	0.2

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85

#### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

25/08/2024 18:27:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)