

CONTAINER FLOVAC 1L (REUSABLE) FOR ASPIRET, ASKIR 20 AND 30



Ref. BCASPI102



25/08/2024 18:28:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CONTAINER FLOVAC 1L (REUSABLE) FOR ASPIRET, ASKIR 20 AND 30



Ref. **BCASPI102**

Description

- FLOVAC® rigid polycarbonate jar / container, reusable.
- Intended to hold FLOVAC disposable pouches- capacity: 1000ml FLOVAC®

® systems therefore replace the use of autoclavable jars, which reduces cleaning and sterilization costs.

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3661040160849	HS Code	4202121900
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 18:28:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr