

CARRYING BAG FOR ASPIRET, ASKIR 20 & ASKIR 30



Ref. BCASPI129



25/08/2024 20:39:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CARRYING BAG FOR ASPIRET, ASKIR 20 & ASKIR 30



Ref. BCASPI129

Description

Suitable for the models below:- New Askir 30 / 30 Proximity- New Askir

230/12V BR- Askir
36 BR- New

Aspiret- New Askir 118 BASIC WM

Features

CE Marking	NONE	Medical device	No
Bar code	3661040165769	HS Code	4202121900
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	21
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	37
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	45
Packed unit / Weight	0.65		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	1	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 20:39:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr